



Recombinant Rat Tubulin polyglutamylase complex subunit 2 (Tpgs2)

Product Code	CSB-MP715258RA
Abbreviation	Tpgs2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AXS8
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MEETPPPMQA GSKPHLEKLT LGVTRILESS PGVTEVSIIE KLPAERHMIS SWEQKNNCVM PEDVRNFYLM TNGFHMTWSV KLDDHIPLG SMVINGISKL TQLIQSSVYS LPNAPTLADL EDDSQEGNED HQLEKPHFDC RSAIFELDSC GGNGKVCLVY KNGKPGLAHD TEIWFLDRAL YWHFLTDTFT AYYRLLITHL GLPQWQYAFT SYGISPQAKQ WFSMYKPITY NTSLLTEEPD TFVNKLDPSK VFKSKNKILI PKKKGPVQPV SGQKGPGLA PPTSKPSAGC GNPVRK
Source	Mammalian cell
Target Names	Tpgs2
Protein Names	Recommended name: Tubulin polyglutamylase complex subunit 2 Short name= PGs2
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.