



Recombinant Rat Protein FAM161A (Fam161a)

Product Code	CSB-EP715264RA-B
Abbreviation	Fam161a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AY14
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MVKLEKMYQD KLTMKDIQAA LVGDDSSSSA SEKSCSHPAL SATSLSEPD DRSSSLSTTT DELPDLEKKT PGEIGTRSYA KELINNMWND FSVEDYTQYD SDLQTAKKNR KKPRAWTPRI TVPVPFEMTV REQKRREKAS DAQETREKML KRNEDDAECK KKFRANPVPS RLLLPLYEDL VKQNEERRKK TRERSKAALL ASQKPFKFA REEQKQAIRE KKLRELCRAK KKPQKFKARP VPRFIYRPPA NVKPKREELY GDSRTQPKAR DVLQSSPWPS HSTYRAFRDP RSPAMPRGKH RHRRLSPSDQ GLEKWKEPFS EQSFRNCPVL CDQCCLYESL CDSNKRQKIL ADIRMGEEIL KETRRPNPSP RHKSPRRSAH ASARPCYSP PMPTASSRGR EQAIRRSEKA RMKELARIGG AR
Source	E.coli
Target Names	Fam161a
Protein Names	Recommended name: Protein FAM161A
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.