



Recombinant Rat ATP-dependent DNA helicase Q1 (Recql)

Product Code	CSB-YP715274RA
Abbreviation	Recql
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AYJ1
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MASIPALTDE LESVSSELHA VDIQIQELTE RQHELLQRKS VLTKRKIQCL EDSAAEASGD CDTSPAAWSK EDFPWSGKVK HVLRDVFKLQ KFRPLQLETV NATMARKDIF LVMPPTGGGKS LCYQLPALCS DGFTLVICPL ISLMEDQLMV LQQLGISATM LNSSSSKEHV KCVHTEMMNK NSHLKLIYVT PEKIAKSKMF MSRLEKAYEA GRLTGAVAVDE VHCCSQWGHG FRPDYKALGI LKRQFPNISL IGLTATATNH VLKDAQKILC VEKCLTFTAS FNRPNLYYEV RQKPSSAEDF IENIANLING RYKGKSGIY CFSQKDSEQV TISLQKLGVR AGTYHANMEP EDRTKVHTQW SANELQVVVA TVAFGMGIDK PDVRFVIHHS MSKSMENYYQ ESGRAGRDDW RADCILYYGF GDIFRISSMV VMENVGQQKL YEMVSYCQNI SKCRRALIAQ HFDEVWNADA CNKMCDNCCK DDSFEKKNIT EHCQALIKIL KQAEGLNEKL TPLKLIDAWM GKGAAKFRVA GVAVPALPRE DLEKIIVHAL LQQYLKEDYS FTAYATISYL KVGPRASLLS NEGHAVTMQV KRSTQSSVRA ASPEACEVDS KGKEKSSAVL C
Source	Yeast
Target Names	Recql
Protein Names	Recommended name: ATP-dependent DNA helicase Q1 EC= 3.6.4.12 Alternative name(s): DNA-dependent ATPase Q1 RecQ protein-like 1
Expression Region	1-621
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the RecQ DNA helicase family. DNA helicases are enzymes involved in various types of DNA repair, including mismatch repair, nucleotide excision repair and direct repair. Some members of this family are associated with genetic disorders with predisposition to malignancy and chromosomal instability. The biological function of this helicase has not yet been



determined. Two alternatively spliced transcripts, which encode the same isoform but differ in their 5 and 3 UTRs, have been described.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.