



Recombinant Rat Secernin-2 (Scrn2)

Product Code	CSB-MP715280RA
Abbreviation	Scrn2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AYR8
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MASSSPDAPC SCDCFVSVPP ASAIPAVIFA KNSDRPRDEV QEVVFIPAGT HTPGSRLQCT YIEVEQVWKT HAVILSRPSW LWGAEMGANE LGVCIGNEAV WTKEPVGEGE ALLGMDLLRL ALERSSTAQA AVHVIAGLLD RYGQGGSCRE DPEPFCYHNT FLLADRTEAW VLETAGRLWA AQRIQGGARN ISNQLSIGTD ISAEHPELRS HAKAQGWWSG QGVFDFAEVF SLTQQPVRME AAKARFRAGC EMLQRQQGSI TAEVMMDILR NKESGICMDS GGFRTTASMV SVLPQDPTKP CVHFLTATPD PSRSVFKPFI FEVGVSQPPQ VLSPTFGAQD PVRILPRFQT QVDRRHSLYR GHQAALGLMD EQEQAQLRQE QQSLEQEGLE ALRGLLTGEQ TPAPQELGSL FQAFVEREHQ VYA
Source	Mammalian cell
Target Names	Scrn2
Protein Names	Recommended name: Secernin-2
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.