



Recombinant Rat Quinone oxidoreductase (Cryz)

Product Code	CSB-MP715281RA
Abbreviation	Cryz
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AYT0
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	ATGQKLMRA IRVFEFGGPE VLKLQSDVVV PAPQSHQVLI K VHACGVNPV ETYIRSGTYS RKPALPYTPG SDVAGIIESV GDGVSAFKKG DRVFCFSTVS GGYAEFALSA DNTTYPLPET LDFRQGAALG IPYFTACRAL FHSARARAGE SVLVHGASGG VGLATCQIAR AHGLKVLGTA GSEEGKKLVL QNGAHEVFNH KEANYIDKIK TSAGDKGVDV IIEMLANKNL SNDLKLLSCG GRVIVVGCRCG SIEINPRDTM AKETSIIGVS LFSSTKEEFQ QFAGILQAGI EKGWVKPVIG SEYPLEKAAQ AHEDIHSSG KMGKMILL
Source	Mammalian cell
Target Names	Cryz
Protein Names	Recommended name: Quinone oxidoreductase EC= 1.6.5.5 Alternative name(s): NADPH:quinone reductase Zeta-crystallin
Expression Region	2-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Crystallins are separated into two classes: taxon-specific, or enzyme, and ubiquitous. The latter class constitutes the major proteins of vertebrate eye lens and maintains the transparency and refractive index of the lens. The former class is also called phylogenetically-restricted crystallins. This gene encodes a taxon-specific crystallin protein which has NADPH-dependent quinone reductase activity distinct from other known quinone reductases. It lacks alcohol dehydrogenase activity although by similarity it is considered a member of the zinc-containing alcohol dehydrogenase family. Unlike other mammalian species, in humans, lens expression is low. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. One pseudogene is known to exist.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.