



Recombinant Rat Mortality factor 4-like protein 1 (Morf4l1)

Product Code	CSB-MP715282RA
Abbreviation	Morf4l1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6AYU1
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MAPKQDPKPK FQEGERVLCF HGPLLYEAKC VKVAIKDKQV KYFIHYSGWN KNWDEWVPES RVLKYVDANL QKQRELQKAN QEQAEGKMR GAAPGKKTSG LQQKNVEVKT KKNKQKTPGN GDGGSTSETP QPPRKKRARV DPTVENEETF MNRVEVKVKI PEELKPWLVD DWDLITRQKQ LFYLPAAKNV DSILEDYANY KKSRTGNTDNK EYAVNEVVAG IKEYFNVMLG TQLLYKFERP QYAEILADHP DAPMSQVYGA PHLLRLFVRI GAMLAYTPLD EKSLALLLNY LHDFLKYLAK NSATLFSASD YEVAPEYHR KAV
Source	Mammalian cell
Target Names	Morf4l1
Protein Names	Recommended name: Mortality factor 4-like protein 1 Alternative name(s): MORF-related gene 15 protein Transcription factor-like protein MRG15
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.