



Recombinant Debaryomyces hansenii DNA replication complex GINS protein PSF1 (PSF1)

Product Code	CSB-EP715329DIS
Abbreviation	PSF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BK12
Product Type	Recombinant Protein
Immunogen Species	Debaryomyces hansenii (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (Torulaspora hansenii)
Purity	≥85% (SDS-PAGE)
Sequence	MYGNLANKLI LDAKRSSNLS ETPLYQTDLV KSIKETTDL NKDAEYLGEE QRIQEEDSEN ADEEKMKINQ CQLFVTHLCM RRNKRCLLAY EKLRAKDIDE FSWLNIDPIT ETGNTDPNKN KFGSTGNGSN SLVSNSANYT TQLNLDNLNH AEQEYYKNYQ HLITNYKSNF ADIDLSGDL D PPTNIFTDVR VLKNGGEVQT EYGVFNLIKD SQFYVRKSDV ERLIQGGYLE EI
Source	E.coli
Target Names	PSF1
Protein Names	Recommended name: DNA replication complex GINS protein PSF1
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.