



Recombinant Debaryomyces hansenii 13 kDa ribonucleoprotein-associated protein (SNU13)

Product Code	CSB-MP715335DIS
Abbreviation	SNU13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BLQ3
Product Type	Recombinant Protein
Immunogen Species	Debaryomyces hansenii (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (Torulaspora hansenii)
Purity	≥85% (SDS-PAGE)
Sequence	MSAPNPKAFP LADSALTQQI LDVVQQSQNL RQLKKGANEA TKTLNRGISE FIIMAADTEP IEILLHLPLL CEDKNVPYVF VPSKTALGRA CGVSRPVIAA SVTTNEASAM KNQIYGKDK IETLLI
Source	Mammalian cell
Target Names	SNU13
Protein Names	Recommended name: 13 kDa ribonucleoprotein-associated protein
Expression Region	1-126
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.