



Recombinant *Debaryomyces hansenii* Pre-mRNA-splicing factor CWC2 (CWC2)

Product Code	CSB-EP715336DIS-B
Abbreviation	CWC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BLU8
Product Type	Recombinant Protein
Immunogen Species	<i>Debaryomyces hansenii</i> (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (<i>Torulasporea hansenii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MSYREPLGEY DCHIDNTRNM PPDTTIVSKS KKGKPARLQV DPESIPDDDR PPQTGNVFNI WFLKWSGGDS STKNYTKSKF RVNIKKDSGY TKAPSNAPLC LFFARGCCYL GKKCSYYHRL PSDTDYFIPT QDCFGRDKTS DYKDDMNGVG SFSKSNRTLY IGGLHMDDKM ENTLTKHFQE FGSIDKIRVL HSKACAFVTF RTENEAQFAK EAMQNSQLDG NEVLNIRWAN EDPNPEAQRQ EKRRLEEVTV NTVKNLLDSV SQTERKTKKV TVEVPDEIEE TESSSEIKAL PSSETSSGLF NNSLNALKQ FQSKKRKIDK PQPKENLPTM LGYSSSDEE
Source	E.coli
Target Names	CWC2
Protein Names	Recommended name: Pre-mRNA-splicing factor CWC2
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.