



Recombinant *Debaryomyces hansenii* Pre-mRNA-splicing factor CWC26 (CWC26)

Product Code	CSB-EP715367DIS
Abbreviation	CWC26
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6BU18
Product Type	Recombinant Protein
Immunogen Species	<i>Debaryomyces hansenii</i> (strain ATCC 36239 / CBS 767 / JCM 1990 / NBRC 0083 / IGC 2968) (Yeast) (<i>Torulasporea hansenii</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MHQLSIFSSL IIGVRMSRAD YLSKYLKDD KKKKKKSKDK KPTNVVVEKN LLDMGDFNDE AEETPVPVIL GAGPKEYKGF KRIDNGDQVS PETQLVTAQD SIPQQETVYR DKSGNIVDIE TKRKQMHEK EQEQIRKRQL EQSINQGDLQ KIEDAEVAAK LDQNKSFVS TRDSEYNELM KSKQRFDDPL QSFSSQKQTV ATSKTGKMYN NKGVSPQNRN NIRAGFFWDG IDRSNGFEDL IVRKRNEVSY TKREASETYD LDMDDDL
Source	<i>E.coli</i>
Target Names	CWC26
Protein Names	Recommended name: Pre-mRNA-splicing factor CWC26
Expression Region	1-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.