



Recombinant *Yarrowia lipolytica* Histone H2A.Z (HTZ1)

Product Code	CSB-MP715402YAP
Abbreviation	HTZ1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6C341
Product Type	Recombinant Protein
Immunogen Species	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MSGGKGKGGK GKGGKGKGGK ASAAVEGAQQ SHSSRAGLQF PVGRIKRYMK KSVQNKVRVG AKAAIYIAAV LEYLTAEVLE LAGNAARDLK VKRITPRHLQ LAIRGDEELD TLIQATIAYG GVMPHINKAL LLKVEQHKKK
Source	Mammalian cell
Target Names	HTZ1
Protein Names	Recommended name: Histone H2A.Z
Expression Region	1-140
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.