



Recombinant *Yarrowia lipolytica* Inositol-pentakisphosphate 2-kinase (IPK1)

Product Code	CSB-EP715415YAP-B
Abbreviation	IPK1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6C6Q7
Product Type	Recombinant Protein
Immunogen Species	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MSSLPTPLPP YKWTLTLAGN ANVVYKSDET DLLLRLRRNR NAPSTAIEVDE YLTGTIRPAI GPFLFHYTVV NLPLGFLESL PEAENLDLGE PLGLLMENLG PKPNETNVLK SHAVKINYSN NWESYTVELK PKWLLQSPTA PKDSINCRTC ALQLKREKPR ICPLKLFNED EQTSLQALD VFPQTQKQFE PLAKFFSNSE LFAEIRHMQH GDELGILGYA NYVQVPPQFV TAMTMRDVSL FVHVQGDSVN GKIVDADLKS VSEKRDYVAS LETDLIEGGW YEKPGTNCLL RN
Source	<i>E.coli</i>
Target Names	IPK1
Protein Names	Recommended name: Inositol-pentakisphosphate 2-kinase EC= 2.7.1.158 Alternative name(s): Inositol-1,3,4,5,6-pentakisphosphate 2-kinase Ins(1,3,4,5,6)P5 2-kinase Short name= InsP5 2-kinase
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.