



Recombinant *Yarrowia lipolytica* RNA polymerase II transcription factor B subunit 3 (TFB3)

Product Code	CSB-BP715419YAP
Abbreviation	TFB3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6C7D0
Product Type	Recombinant Protein
Immunogen Species	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>)
Purity	≥85% (SDS-PAGE)
Sequence	MSEDGDICPI CKSSRYLNPD MKFLVNPQCY HKMCESCVDL LFAYGPVTCP HNGCEKILRK NKFQKQIFED VAVEKEVDVR QRVMSVMNKR EDEFDTLNDY NAYLEKIEDS IFTLLNGTAD EKDALTKIEI AYEKEHKAEI IANNKLRDEE DKYEQKEDW MKEQRKANYL MAVEEQRAAQ EDKEAAKREL VHQLATSNKD ASSIEQNVQK TALKRSSART VNFQPMPTSR YFEKVKKESS KPETPFTPFN GDRQIPALFE VQDHYDDKVI EPLAHDIQHQ AGGFTLKAAY ARILQQAFFG LGVDVQTELL AEENGNSNGT AVPVKTEEDP VKLEEPSIQV KVESTV
Source	Baculovirus
Target Names	TFB3
Protein Names	Recommended name: RNA polymerase II transcription factor B subunit 3 Alternative name(s): RNA polymerase II transcription factor B 38 kDa subunit RNA polymerase II transcription factor B p38 subunit
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.