



Recombinant *Yarrowia lipolytica* Phosphatidylglycerol/phosphatidylinositol transfer protein (NPC2)

Product Code	CSB-BP715454YAP
Abbreviation	NPC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CHU5
Product Type	Recombinant Protein
Immunogen Species	<i>Yarrowia lipolytica</i> (strain CLIB 122 / E 150) (Yeast) (<i>Candida lipolytica</i>)
Purity	>85% (SDS-PAGE)
Sequence	DEDIPGDS PISLCDAEYE QLLEIKHLSI DPNPPAKGQN LTIEASGYLY EDVEEGAYIE VEVRYGYIRL ISQTLDLCEQ SEQVDWSCPI KAGDLKLNKQ IELPNEIPPG KYVAVARAYT VDDDLITCLL TTVTFPASL
Source	Baculovirus
Target Names	NPC2
Protein Names	Recommended name: Phosphatidylglycerol/phosphatidylinositol transfer protein Short name= PG/PI-TP
Expression Region	53-189
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.