



Recombinant Kluyveromyces lactis Protein IBD2 (IBD2)

Product Code	CSB-BP715505KBK
Abbreviation	IBD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CRM3
Product Type	Recombinant Protein
Immunogen Species	Kluyveromyces lactis (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (Candida sphaerica)
Purity	≥85% (SDS-PAGE)
Sequence	MGNGNASIEF VSENGPIDFS VMMQEGVKAL TKILSTHLQD NPDLINEKSM KVIKDNNGK LEPVFAPSEG NSGLINGNDQ DENYDEDRDI EEANEDYEVI EFDNEGNYHV CAKGTQNRAI GYDEKVYSKS EHYADSTDTH DYDEDGIINE DNNQLLPKDK WNGHMTRLKN GAISSAEHEV VFEHDNGDNT TFSPAETIND GHMDKRDRSH MHHSQHSEKT GGSCACQENK DSFKYDPHVP NSAPNFRALL QHTVNGKQMC PFCEYYVVFV TPPAYIMRWA SSKISAEHSQ YKMNDQIVDG DSGDMDRNL D SVETNEMNKQ QTKKKKKKSK RSKR
Source	Baculovirus
Target Names	IBD2
Protein Names	Recommended name: Protein IBD2
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.