



Recombinant Kluyveromyces lactis RNA polymerase II transcription factor B subunit 4 (TFB4)

Product Code	CSB-MP715530K BK
Abbreviation	TFB4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6CVX9
Product Type	Recombinant Protein
Immunogen Species	Kluyveromyces lactis (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) (Candida sphaerica)
Purity	≥85% (SDS-PAGE)
Sequence	MDAIADTVFQ GTKSKTQEIE DTPSLLTVVV DTSIHSWVQL TKQQSGSGSE GSSGEKQLIE ALKSIVVFLN AHLAFNSGNQ VCLIAAHSEG MKYLYPSADS KPSMSMVSSD MYRGFRNVDE IVVEQWYRLF KEELEGQESK VSMKSSLSGA MSSALTYVNR ILKENENTSL RSRLLVITCG TSQGKDEIFQ YIPIMNCIFS ATKMKCSIDV VKIGGGIEST FLQQATDATS GVYLHVENTR GLIQYLSTAM FIDPSLRNVI IKPNQGSVDF RTSCFLTGV VAVGFVCSVC LCVLSVIPPG QKCPACDSPF DSKIIARLRR KPVLSNGV PK KKNVSNK
Source	Mammalian cell
Target Names	TFB4
Protein Names	Recommended name: RNA polymerase II transcription factor B subunit 4 Alternative name(s): RNA polymerase II transcription factor B 34 kDa subunit RNA polymerase II transcription factor B p34 subunit
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.