



Recombinant Erwinia carotovora subsp. atroseptica Ribosomal RNA small subunit methyltransferase H (rsmH)

Product Code	CSB-MP715565EAAB
Abbreviation	rsmH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6D0H5
Product Type	Recombinant Protein
Immunogen Species	Pectobacterium atrosepticum (strain SCRI 1043 / ATCC BAA-672) (Erwinia carotovora subsp. atroseptica)
Purity	≥85% (SDS-PAGE)
Sequence	MLENYKHTTV LLDEAVNGLN IRSGGIYIDG TFGRGGHSRL ILSQLGPEGR LLAIDRDPQA IEAAKAIDDS RFSIIHGPFs AMADYVAELG LTGQIDGVLL DLGVSSPQLD DPERGFSFMR DGPLDMRMDP TRGSSAAEWL MKAEADDIVW VLKTFGEERF AKRIARAIVE RNRTEPMTRT KELASLIAAA SPIREKHKHP ATRSFQAIRI YINSELEEIE RALEGALSVL APQGRLSVIS FHSLEDRIK RFIRHQSRGP QVPAGLPLTE EQLRSQGGQT LKPVGKKLMP SEAEVAENPR ARSSVLRFAE RLPA
Source	Mammalian cell
Target Names	rsmH
Protein Names	Recommended name: Ribosomal RNA small subunit methyltransferase H EC= 2.1.1.199 Alternative name(s): 16S rRNA m(4)C1402 methyltransferase rRNA (cytosine-N(4)-)-methyltransferase RsmH
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.