



Recombinant *Erwinia carotovora* subsp. atroseptica UDP-4-amino-4-deoxy-L-arabinose--oxoglutarate aminotransferase (arnB)

Product Code	CSB-EP715581EAAB
Abbreviation	arnB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6D2E9
Product Type	Recombinant Protein
Immunogen Species	Pectobacterium atrosepticum (strain SCRI 1043 / ATCC BAA-672) (<i>Erwinia carotovora</i> subsp. atroseptica)
Purity	≥85% (SDS-PAGE)
Sequence	MTDFLPFSRP TMGDEEIAAV AEVLRSGWIT TGPKCQQLEQ AFCQQFGCQQ AIAVSSATGG MHVTLMALGI GPGDEVITPS QTWVSTVNII TLLGAEPVMV DVDRHTLMVR SQDIEAAITP KTKAIIPVHY AGAPADLDAL RALSERYGIP LIEDAAHAVG TQYRGEWIGA RGTAIFSFA IKNITCAEGG MVVTDDEALA ERIRSLKFHG LGVDAFDRQL QGRKPQAEVV TPGFKYNLAD INAAIALVQL DKLPAINARR QQLAARYLTQ LHSLPLQPLA VPDYPHLHAW HLFMVRVDET QCGISRDDL M AALQTHGIGT GLHFRAVHTQ KYRERYPHL HLPETEWNSA SLMTLPLFPD MQDSDVDRVV AALTSILESV RD
Source	E.coli
Target Names	arnB
Protein Names	Recommended name: UDP-4-amino-4-deoxy-L-arabinose--oxoglutarate aminotransferase EC= 2.6.1.87 Alternative name(s): UDP-(beta-L-threo-pentapyranosyl-4"-ulose diphosphate) aminotransferase Short name= UDP-Ara4O aminotransferase UD
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.