



Recombinant *Xenopus laevis* Abhydrolase domain-containing protein FAM108B1 (fam108b1)

Product Code	CSB-MP715644XBE
Abbreviation	fam108b1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DCC5
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	KIASKLAFLP PDPTYTLICD ESGSRWTLHL SERADWQYSS REKDAIECFM TRTSRGNRIA CMFVRCCPSA KYTLLEFSHGN AVDLGQMSSF YIGLGSRINC NIFSVDYSGY GSSSGKPSEK NLYADIDAAW IALRTRYGIR PEHVIIYGQS IGTVPSVDLA ARYESAAVIL HSPLTSGMRV AFPDTKKTYC FDAFPNIDKI SKITSPVLII HGTEDEVIDF SHGLALFERC QRPVEPLWVE GAGHNDVELY GQYLERLKQF VAQELVNL
Source	Mammalian cell
Target Names	abhd17b
Protein Names	Recommended name: Abhydrolase domain-containing protein FAM108B1 EC= 3.-.-
Expression Region	21-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.