



Recombinant Danio rerio Pre-mRNA-splicing factor 38A (prpf38a)

Product Code	CSB-MP715687DIL
Abbreviation	prpf38a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6DHU4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MANRTVKDAN SVHGTPNPQYL VEKIIRTRIIY ESKYWKEECF GLTAEVVVDK AMELKFGGV YGGNVKPTPF LCLTLKMLQI QPEKDIIVEF IKNEDFKYVR LLGAMYMRLT GTSVDCYKYL EPLYNDYRKI KSQNRNGEFE LMHVDEFIDE LLHAERMCDI ILPRLQKRQV LEEAELL DPR ISALEEDLDE VETSEEEEDDE DEKPERMQSP EPHRRSYRDM DRPRRSPSPR YRRSRSPRRR SRSPKRRSPS PRRERDRDRH RSKSPRRHRS RSRERRHRSK SPGHRSHRH RSHSKSPESR SKKSHKRSRR GNE
Source	Mammalian cell
Target Names	prpf38a
Protein Names	Recommended name: Pre-mRNA-splicing factor 38A
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.