



# Recombinant Influenza A virus Nucleoprotein (NP)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-MP715711IBM                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Abbreviation</b>      | NP                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                            |
| <b>Uniprot No.</b>       | Q6DPE7                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Immunogen Species</b> | Influenza A virus (strain A/Chicken/Hong Kong/37.4/2002 H5N1 genotype X2)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Sequence</b>          | MASQGTKRSY EQMETGGERQ NATEIRASVG RMVSGIGRFY IQMCTELKLS<br>DYEGRLIQNS ITIERMVLSA FDERRNRYLE EHPSAGKDPK KTGGPIYRRR<br>DGKWVRELIL YDKEEIRRIW RQANNGEDAT AGLTHLMIWH SNLNDATYQR<br>TRALVRTGMD PRMCSLMQGS TLPRRSGAAG AAIKGVGTMV MELIRMIKRG<br>INDRNFWRGE NGRRTRIAYE RMCNILKGKF QTAAQRAMMD QVRESRNPNGN<br>AEIEDLIFLA RSALILRGSV AHKSCLPACV YGLAVASGYD FEREGYSLVG<br>IDPFRLQNS QVFSLIRPNE NPAHKSQLVW MACHSAAFED LRVSSFIRGT<br>RVVPRGQLST RGVQIASNEN MDTMDSNTLE LRSRYWAIRT RSGGNTNQQR<br>ASAGQISVQP TFSVQRNLPF ERATIMAAFT GNTEGRTSDM RTEIIRMMES<br>ARPEDLSFQG RGVFELSDEK ATNPVPSFD MSNEGSYFFG DNAEEFES |
| <b>Source</b>            | Mammalian cell                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Target Names</b>      | NP                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Protein Names</b>     | Recommended name: Nucleoprotein Alternative name(s): Nucleocapsid protein<br>Short name= Protein N                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Expression Region</b> | 1-498                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                            |