



# Recombinant *Oryza sativa* subsp. japonica B3 domain-containing protein Os02g0683500 (Os02g0683500, LOC\_Os02g45850)

<b>Product Code</b>	CSB-EP715761OFG-B
<b>Abbreviation</b>	Os02g0683500, LOC_Os02g45850
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6EU30
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEFTTSSRFS KEEEDDEEQDE AGRREIPFMT ATAEAAPAPT SSSSSPAHHA ASASASASAS GSSTPFRSDD GAGASGSGGG GGGGGAEV EKEHMFDPKVV TPSDVGKLNLR LVIPKQYAEK YFPLDAAANE KLLLLNFEDR AGKPWRFRYS YWNSSQSYVM TKGWSRFVKE KRLDAGD TVS FSRGIGDEAA RHRLFIDWKR RADTRDPLRL PRGLPLP MPL TSHYAPWGIG GGGGFFVQPS PPATLYEHRL RQGLDFRAFN PAAAMGRQVL LFGSARIPPQ APLLARAPSP LHHHYTLQPS GDGVRAAGSP VVLD SVPVIE SPTTAAKRVR LFGVNLDNPH AGGGGGAAAG ESSNHGNALS LQTPAWMRRD PTLRLLLELPP HHHHGAESSA ASSPSSSSSS KRDAHSALDL DL
<b>Source</b>	E.coli
<b>Target Names</b>	Os02g0683500
<b>Protein Names</b>	Recommended name: B3 domain-containing protein Os02g0683500
<b>Expression Region</b>	1-412
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.