



Recombinant Salmonella paratyphi A ADP-L-glycero-D-manno-heptose-6-epimerase (hldD)

Product Code	CSB-YP715820SBM
Abbreviation	hldD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PC05
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	MIIVTGGAGF IGSNIVKALN DKGITDILVV DNLKDGTKFV NLVDLNIADY MDKEDFLIQI MSGEELGDIE AIFHEGACSS TTEWDGKYMM DNNYQYSKEL LHYCLEREIP FLYASSAATY GGRTSDFIES REYEKPLNVY GYSKFLFDEY VRQILPEANS QIVGFYFNV YGPREGHKGS MASVAFHLNT QLNNGESPKL FEGSENFKRD FVYVGDVAAV NLWFLESGKS GIFNLGTGRA ESFQAVADAT LAYHKKGSIE YIPFPDKLKG RYQAFQTADL TNLRNAGYDK PFKTVAEGVT EYMAWLNLRDA
Source	Yeast
Target Names	hldD
Protein Names	Recommended name: ADP-L-glycero-D-manno-heptose-6-epimerase EC=5.1.3.20 Alternative name(s): ADP-L-glycero-beta-D-manno-heptose-6-epimerase Short name= ADP-glyceromanno-heptose 6-epimerase Short name= ADP-hep 6-epimerase
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.