



Recombinant Salmonella paratyphi A FMN-dependent NADH-azoreductase (azoR)

Product Code	CSB-MP715870SBM
Abbreviation	azoR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PHS8
Product Type	Recombinant Protein
Immunogen Species	Salmonella paratyphi A (strain ATCC 9150 / SARB42)
Purity	≥85% (SDS-PAGE)
Sequence	MSKVLVLKSS ILAGYSQSGQ LTDYFIEQWR EKHVADEITV RDLAANPVPV LDGELVGAMR PGDAPLTPRQ QDALALSDEL IAEKHAHDVI VIAAPMYNFN IPTQLKNYFD LIARAGITFR YTEKGPEGLV TGKRAVLSS RGGIHKDTP DLIAPYLKVF LGFIGITDVN FVFAEGIAYG PEVAAKAQAD AKAAIDSVVA A
Source	Mammalian cell
Target Names	azoR
Protein Names	Recommended name: FMN-dependent NADH-azoreductase EC= 1.7.-.- Alternative name(s): Azo-dye reductase FMN-dependent NADH-azo compound oxidoreductase
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.