



Recombinant *Xenopus laevis* Eukaryotic translation initiation factor 3 subunit H (eif3h)

Product Code	CSB-MP715936XBE
Abbreviation	eif3h
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5PPY6
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	≥85% (SDS-PAGE)
Sequence	MAARKETGNT APIAETAVKQ VQIDGLVVLK IIKHYQEEGH GSEVVQGVLL GLVVDDRLEI TNCFPFPQHT EDDVDFDEVQ YQMEMMRSLR HVNIDHLHVG WYQSTFYGTF VSRALLDSQF SYQHAIEESV VLIYDPIKTS QGSLSLKAYR LTPKLMEICK EKDFS AGLK KANVAYEDMF EEVPIVIKNS HLINVLVWEL DKKAPVTEKH ELLNLSSSNH LEKSLQLLMD RVDEMSQDIV KYNTYQRNVG KQQQKQKHQYT QRKQQENLQR LSRGETPLPE EDVNKMFKPP VTPSRMDPLL IAGQINTYIQ HIKGFTSQNL GKL FMAEALQ DQTN
Source	Mammalian cell
Target Names	eif3h
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit H Short name= eIF3h Alternative name(s): Eukaryotic translation initiation factor 3 subunit 3 eIF-3-gamma eIF3 p40 subunit
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.