



# Recombinant Fundulus heteroclitus Serine/threonine-protein kinase Sgk1 (sgk1)

<b>Product Code</b>	CSB-YP715951FDK
<b>Abbreviation</b>	sgk1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5Q0U5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Fundulus heteroclitus (Killifish) (Mummichog)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTIKTETEKP ALTYSKTRGL VELITAFMKQ RRMGLNDFIQ KLATNSYACK HPEVQSILNL TPPQDPELMN SNPSPPSPS QQINLGPSSN PSAKPSDFHF LKVIGKGSFG KVLARHRTD DQFYAVKVLQ KKAILKKKEE KHIMSERNVL LKNVKHPFLV GLHYSFQTAD KLYFVLDYIN GGELFYHLQR ERCFLEPRAR FYSAEIASAL GYLHSLNIVY RDLKPENILL DSQGHILTD FGLCKENIEP NGTTSTFCGT PEYLAPEVLH KQPYDRTVDW WCLGAVLYEM LYGLPPFYSR NTAEMYDNIL NKPLQLKPNL SNAARHLEGG LLQKDRTRKRL GCKDDFTEIK NHVFFSPINW DDLNAKKMTP PFNPNVTGPN DLRHFDPEFT DEPVPSSIGC SPDCALATAS IKEAAEAFVG FSYAPSMDSY L
<b>Source</b>	Yeast
<b>Target Names</b>	sgk1
<b>Protein Names</b>	Recommended name: Serine/threonine-protein kinase Sgk1 EC= 2.7.11.1 Alternative name(s): Serum/glucocorticoid-regulated kinase 1
<b>Expression Region</b>	1-431
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.