



Recombinant Pongo abelii Nuclear receptor-binding factor 2 (NRBF2)

Product Code	CSB-MP716027PYX
Abbreviation	NRBF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R4C9
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MEVMEGPLNL AHQQSRRADR LLAAGKYEEA ISCHKKAAAY LSEAMKLTQS EQAHLSELELQ RDSHMKQLLL IQRWKRAQR EERLKAQQNT DKDVAAHLQT SHKPSAEDAE GQSPLSQKYS PSTEKCLPEI QGIFDRDPDT LLYLLQQKSE PAEPCIGSKA PKDDKTIIIE QATKIADLKR HVEFLVAENE RLRKENKQLK AEKARLLKGP IEKELDVAD FVETSELWSL PPHSETATAS STWQKFAANT GKAKDIPIP LPPLDFPSPE LPLMELSEDI LKGFMMN
Source	Mammalian cell
Target Names	NRBF2
Protein Names	Recommended name: Nuclear receptor-binding factor 2 Short name= NRBF-2
Expression Region	1-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.