



Recombinant Pongo abelii Ankyrin repeat and SOCS box protein 6 (ASB6)

Product Code	CSB-MP716032PYX
Abbreviation	ASB6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R4M7
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MPFLHGFRRRI IFEYQPLVDA ILGSLGIQDP ERQESLDRPS YVASEESRIL VLTELLERKA HSPFYQEGVS NALLKMAELG LTRAADVLLR HGANLNFEDP VTYYTALHIA VLRNQPDME LLVHHGADIN RRDRIHESSP LDLASEEPER LPCLQRLLDL GADVNAADKH GKTALLHALA SSDGVQIHNT ENIRLLLEGG ADVKATTKDG DTVFTCIIFL LGETVGGDKE EAQMINRFCF QVTRLLLAHG ADPSECPAHE SLTHICLKSF KLHFPLLRFL LESGAAYNCS LHGASCWSGF HIIFERLCSH PGCTEDESHT DLLRKAETVL DLMVTNSQKL QLPENFDIHP VGSLAEKIQA LHFSLRQLES YPPPLKHLCR VAIRLYLQPW PVDVKVKALP LPDRWKWYLL SEHSGSMEDD I
Source	Mammalian cell
Target Names	ASB6
Protein Names	Recommended name: Ankyrin repeat and SOCS box protein 6 Short name= ASB-6
Expression Region	1-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.