



Recombinant Pongo abelii Alpha-1-antichymotrypsin (SERPINA3)

Product Code	CSB-BP716041PYX
Abbreviation	SERPINA3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R536
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	HPNCPLD EENPTQENQD RGTHVDLGLA STNVDFAFSL YKQLVLKAPD KNVIFSPLSI STALAFSLG AHNTTLTEIL TGLRFNLTET SEAEIHQSFQ HLLRTLNQSS DELQLSMGNA MFVEEQLSLL DRFMEDAKRL YGSEAFATDF QDSAAAKKLI NDYVKNRTRG KITDLIKDLD SQTMMVLVNY IFFKAKWKMP FDPQDTHQSR FYLSKKKWVM VPMMSLHHLT TPYFRDEELS CTVVELKYTG NASALFILPD QDKMEEVEAM LLPETLKRWR DSLEFRRIDE LYLPKFISISR AFNLENILLQ LGIVEAFTSK ADLSGITGAR NLVVSQVVHK AVLDVFEEGT EASAATAVKI TLLSALVDPM TIVRFNRPFL MIIVPTDTQN LLFISKVINP KQA
Source	Baculovirus
Target Names	SERPINA3
Protein Names	Recommended name: Alpha-1-antichymotrypsin Short name= ACT Alternative name(s): Serpin A3
Expression Region	24-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.