



Recombinant Pongo abelii GTP-binding protein Di-Ras2 (DIRAS2)

Product Code	CSB-EP716069PYX
Abbreviation	DIRAS2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5R6S2
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MPEQSNDYRV AVFGAGGVGK SSLVLRVFKG TFRESYIPTV EDTYRRVISC DKSICTLQIT DTTGSHQFPA MQRLSISKGH AFILVYSITS RQSLEELKPI YEQICEIKGD VESIPIMLVG NKCDSPSRE VQSSEAEALA RTWKCAFMET SAKLNHNVKE LFQELLNLEK RRTVSLQIDG KKSQQQRKE KLKGKC
Source	E.coli
Target Names	DIRAS2
Protein Names	Recommended name: GTP-binding protein Di-Ras2 Alternative name(s): Distinct subgroup of the Ras family member 2
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.