



Recombinant Pongo abelii Perilipin-3 (PLIN3)

Product Code	CSB-EP716146PYX
Abbreviation	PLIN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RAV8
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	SADGAEADG STQVTVEEPV QQPSVVDRA SMLPISSTCD MVSAAYASTK ESYPHVKTVC DAAEKGVRTL TAAAVSGAQP ILSKLEPQIA SASEYHRGL DKLEENLPIL QQPTERVLAD TKELVSSKVS GAQEMVSSAK DTVATQLSEA VDATRGAVQS GVDKTKSVVT GGVQSVMGSR LGQMVLSGVD TVLGKSGEWA DNHLPLTDAE LARIATSLDG FDVASVQQQR QEQSYFVRLG SLSERLRQHA YEHS LGKLRA TKQRAQEALL QLSQALSLME TVKEGVDQKL VEGQEKLHQM WLSWNQKQLQ GPEKEPPKPE QVESRALTMF RDIAQQLQAT CTSLGSSIQG LPTNVKDQVQ QARRQVEDLQ ATFSSIHSFQ DLSSSILAQS RERVASAREA LDHMVEYVAQ NTPVTWLVGP FAPGITEKAP EEKK
Source	E.coli
Target Names	PLIN3
Protein Names	Recommended name: Perilipin-3 Alternative name(s): Mannose-6-phosphate receptor-binding protein 1
Expression Region	2-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.