



Recombinant Pongo abelii Protein NDRG2 (NDRG2)

Product Code	CSB-EP716160PYX
Abbreviation	NDRG2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RBN6
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	AELQEVQIT EEKPLLPGQT PEAAKEAELA ARILLDQGQT HSVETPYGSV TFTVYGTPKP KRPAITYHD VGLNYKSCFQ PLFQFEDMQE IIQNFVRVHV DAPGMEEGAP VFPLGYQYPS LDQLADMIPC VLQYLNFTSI IGVGVGAGAY ILARYALNHP DTVEGLVLIN IDPNAKGWMD WAAHKLTGLT SSIPEMILGH LFSQEELSGN SELIQKYRNI ITHAPNLNI ELYWNSYNNR RDLNFERGGD ITLKCPVMLV VGDQAPHEDA VVECNSKLDP TQTSFLKMAD SGGQPQLTQP GKLTEAFKYF LQGMGYMASS CMTRLSRSRT ASLTAASVD GNRSRRTLS QSSESGTLSS GPPGHTMEVS C
Source	E.coli
Target Names	NDRG2
Protein Names	Recommended name: Protein NDRG2 Alternative name(s): N-myc downstream-regulated gene 2 protein
Expression Region	2-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.