



Recombinant Pongo abelii Properdin (CFP)

Product Code	CSB-EP716161PYX-B
Abbreviation	CFP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RBP8
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	DPV LCFTRYEESS GKCKGLLGGG VSVKDCCLNT ASAYQERSGG LCQPCRSRW SLWSTWAPCS VTCSEGSQRLR YRRCVGNWNGQ CSEKVALGTL EWQLQACEDQ QCCPEMGGWS GWGPWEPCSV TCSKGTRTRR RACNHPAPKC GGHCPGQAQE SEACDTQQSC PTHGAWAAWG PWTPCSGSCH SGTHEPKETR SRKCSAPEPS QKPPGKPCPG LAYEQRRCTG LPPCPVAGSW GPWGPVSPCP VTCGLGQTME QRTCDHPVPQ HGGPFCAGDA TRTHICNTAV PCPVDGEWDS WGEWSTCVRR NMKSISCQEI PGQQSRWRTC KGRKFDGHRC AGQQQDIRHC YSIQHCPKLG SWSEWSTWGL CIPPCGPNPT RARQRLCTPL LPKYPPTVSM VEGQGEKNVT FWGRPLPRCE ELQGQKLVVE EKRPCLHVPA CKDPEEEKL
Source	E.coli
Target Names	CFP
Protein Names	Recommended name: Properdin Alternative name(s): Complement factor P
Expression Region	28-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.