



Recombinant Pongo abelii DNA replication complex GINS protein PSF3 (GINS3)

Product Code	CSB-BP716205PYX
Abbreviation	GINS3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RDV0
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	MSEAYFRVES GALGPEENFL SLDDILMSHE KLPVRTETAM PRLGAFFLER SAGAETDNAV PQGSKLELPL WLAKGLFDNK RRILSVELPK IYQEGWRTVF SADPNVVDLH KMGPHFYGFG SQLLHFDSPE NADISQSLQ TFIGRFRRIM DSSQNAYNED TSALVARLDE MERGLFQIGQ NGLNDFQCWE KGQASQITAS NLVQNYKKRK FTDMED
Source	Baculovirus
Target Names	GINS3
Protein Names	Recommended name: DNA replication complex GINS protein PSF3 Alternative name(s): GINS complex subunit 3
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.