



# Recombinant Pongo abelii Selenium-binding protein 1 (SELENBP1)

<b>Product Code</b>	CSB-BP716233PYX
<b>Abbreviation</b>	SELENBP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5RF48
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ATKCGNCGP GYSTPLEAMK GPREEIVYLP CIYRNTGTEA PDYLATVDVD PKSPQYCQVI HRLPMPNLKD ELHHSWNTC SSCFGDSTKS RTKLVLP SLI SSRIYVVDVG SEPRALKLHK VIEPKDIHAK CELAFLHTSH CLASGEVMIS SLGDVKGNGK GGFVLLDAET FEVKGTWERR GGAAPLG YDF WYQPRHNMI STEWAAPNVL RDGFNPADVE AGLYGSHLYV WDWQRHEIVQ TSLKDG LIP LEIRFLHNPD AAQGFVGCAL SSTIQRFYKN EGGTWSVEKV IQVPPKVKKG WLLPEMPGLI TDILLSLDDR FLYFSNWLHG DLRQYDISDP QRPRLTGQLF VGG SIVKGGP VQVLEDQELK SQPEPLVVKG KRVAGGPQMI QLSLDGKRLY VTTSLYSAWD KQFY PD LIRE GSVMLQVDVD TVKGGLK LNP NFLVDFGKEP LGPALAHEL R YPGGDCSSDI WI
<b>Source</b>	Baculovirus
<b>Target Names</b>	SELENBP1
<b>Protein Names</b>	Recommended name: Selenium-binding protein 1
<b>Expression Region</b>	2-472
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.