



Recombinant Pongo abelii Selenium-binding protein 1 (SELENBP1)

Product Code	CSB-EP716233PYX-B
Abbreviation	SELENBP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RF48
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	≥85% (SDS-PAGE)
Sequence	ATKCGNCGP GYSTPLEAMK GPREEIVYLP CIYRNTGTEA PDYLATVDVD PKSPQYCQVI HRLPMPNLKD ELHHSWNTC SSCFGDSTKS RTKLVLP SLI SSRIYVVDVG SEPRALKLHK VIEPKDIHAK CELAFLHTSH CLASGEVMIS SLGDVKGN GK GGFVLLDAET FEVKGTWERR GGAAPLG YDF WYQPRHNMI STEWAAPNVL RDGFNPADVE AGLYGSHLYV WDWQRHEIVQ TSLKDG LIP LEIRFLHNPD AAQGFVGCAL SSTIQR FYKN EGGTWSVEKV IQVPPKVKG WLLPEMPGLI TDILLSLDDR FLYFSNWLHG DLRQYDISDP QRPRLTGQLF VGG SIVKGGP VQVLEDQELK SQPEPLVVKG KRVAGGPQMI QLSLDGKRLY VTTSLYSAWD KQFY PD LIRE GSVMLQVDVD TVKGGLK LNP NFLVDFGKEP LGPALAHEL R YPGGDCSSDI WI
Source	E.coli
Target Names	SELENBP1
Protein Names	Recommended name: Selenium-binding protein 1
Expression Region	2-472
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.