



Recombinant Cyclin-D (CycD)

Product Code	CSB-EP716272OCT
Abbreviation	CycD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SCB5
Product Type	Recombinant Protein
Immunogen Species	Ostreococcus tauri
Purity	≥85% (SDS-PAGE)
Sequence	MNRSRTSSFS TDRSASTVAL APDGDYERER VSALASWPEY AVDALRSEQD FIVHQSAALW TH RTPALHDP DSQETQLDVA KGLLARERET HGSFVFDARA AHHCAFRSQL VEWILDVCAG ERFGPTTADV AIAYTDRVLS KTVVPKTS LH LVALCCLHIA VKYEEIEERV PTMSKLR SWT SNMYSPEIIR KMELAVLIEL GWD LGVLT PA HFLESFLALT NGGISD GDDI EHGDAYKERY REELRYFVCQ LYSLCVQDTS LLNQPPSQIA SAVIATARVH LGVKPMCSPE LRAAGNVTPQ QIYPLVAHML KLWDEACAED EAMDEVETSA EFNSLTIQVP KPIGHDIVSK MGYEHVRTEN ASPTCPFDMQ WDEE
Source	E.coli
Target Names	CycD
Protein Names	Recommended name: Cyclin-D
Expression Region	1-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.