



Recombinant *Thermus thermophilus* 1,4-alpha-glucan branching enzyme TTHA1902 (TTHA1902), partial

Product Code	CSB-YP716285TNT
Abbreviation	TTHA1902
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SH28
Product Type	Recombinant Protein
Immunogen Species	<i>Thermus thermophilus</i> (strain HB8 / ATCC 27634 / DSM 579)
Purity	≥85% (SDS-PAGE)
Source	Yeast
Target Names	TTHA1902
Protein Names	Recommended name: 1,4-alpha-glucan branching enzyme TTHA1902 EC= 2.4.1.18 Alternative name(s): 1,4-alpha-D-glucan:1,4-alpha-D-glucan 6-glucosyl-transferase Alpha-(1->4)-glucan branching enzyme Branching enzyme Short name= BE
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.