



Recombinant Mouse Luc7-like protein 3 (Luc7I3)

Product Code	CSB-YP716334MO
Abbreviation	Luc7I3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SUF2
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MISAAQLLDE LMGRDRNLAP DEKRSNVRWD HESVCKYYLC GFCPAELFTN TRSDLGPCEK IHDENLRKQY EKSSRFMKVG YERDFLYLQ SLLAEVERRI RRGHARLALS QNQQSSGAAG PTGKNEEKIQ VLTDKIDVLL QQIEELGSEG KVEEAQGMMK LVEQLKEERE LLRSTTSTIE SFAAQEKQME VCEVCGAFLI VGDAQSRVDD HLMGKQHMGY AKIKATVEEL KEKLRKRTEE PDRDERLKKE KQEREEREKE REREREERER KRRREEEERE KERARDRERR KRSRSRSRHS SRTSDRRCR SRDHKRSRSR DRRRSRSDR RRRSRSHDRSE RKHRSRSDR RRSKSRDRKS YKHRKSRDR EQDRKSKEKE KKGSDDKKSS VKSSSREKQS EDTNPESKES DTKNEVNGTS EDIKSEGDTQ SN
Source	Yeast
Target Names	Luc7I3
Protein Names	Recommended name: Luc7-like protein 3 Alternative name(s): Cisplatin resistance-associated-overexpressed protein
Expression Region	1-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.