



# Recombinant Human PH domain-containing protein C10orf81 (C10orf81)

<b>Product Code</b>	CSB-BP716341HU
<b>Abbreviation</b>	C10orf81
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5SXH7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAGGKQFTFS YENEVCKQDY FIKSPPSQLF SSVTSWKKRF FILSKAGEKS FSLSYKDDHH HRGSIEIDQN SSVEVGISSQ EKMQSVQKMF KCHPDEVMSI RTTNREYFLI GHDREKIKDW VSFMSSFRQD IKATQQNTEE ELSLGNKRTL FYSSPLLGPS STSEAVGSSS PRNGLQDKHL MEQSSPGFRQ THLQDLSEAT QDVKEENHYL TPRSVLLELD NIIASSDSGE SIETDGPDQV SGRIECHYEP MESYFFKETS HESVDSSKEE PQTLPETQDG DLHLQEQQSG IDWCLSPADV EAQTTNDQKG SASLTVVQLS ILINNIPDES QVEKLNVFLS PPDVINYAL TEATGRICVS QWEGPPRLGC IFCHGDHLLA VNDLKPQSLE EVSLFLTRSI QKEKLTIG RIPNSETFHA ASCMCPSKCQ SAAPSQLDKP RLNRAPKRSP AIKKSQQKGA RE
<b>Source</b>	Baculovirus
<b>Target Names</b>	PLEKHS1
<b>Protein Names</b>	Recommended name: PH domain-containing protein C10orf81 Alternative name(s): Epididymis luminal protein 185 Short name= hEL185
<b>Expression Region</b>	1-462
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.