



Recombinant Human PH domain-containing protein C10orf81 (C10orf81)

Product Code	CSB-EP716341HU
Abbreviation	C10orf81
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5SXH7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAGGKQFTFS YENEVCKQDY FIKSPPSQLF SSVTSWKKRF FILSKAGEKS FSLSYKDDHH HRGSIEIDQN SSVEVGISSQ EKMQSVQKMF KCHPDEVMSI RTTNREYFLI GHDREKIKDW VSFMSSFRQD IKATQQNTEE ELSLGNKRTL FYSSPLLGPS STSEAVGSSS PRNGLQDKHL MEQSSPGFRQ THLQDLSEAT QDVKEENHYL TPRSVLLELD NIIASSDSGE SIETDGPDQV SGRIECHYEP MESYFFKETS HESVDSSKEE PQTLPETQDG DLHLQEQQSG IDWCLSPADV EAQTTNDQKG SASLTVVQLS ILINNIPDES QVEKLNVFLS PPDVINYAL TEATGRICVS QWEGPPRLGC IFCHGDHLLA VNDLKPQSLE EVSLFLTRSI QKEKLTIG RIPNSETFHA ASCMCPSKCQ SAAPSQLDKP RLNRAPKRSP AIKKSQQKGA RE
Source	E.coli
Target Names	PLEKHS1
Protein Names	Recommended name: PH domain-containing protein C10orf81 Alternative name(s): Epididymis luminal protein 185 Short name= hEL185
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.