



Recombinant Human Uncharacterized protein C1orf111 (C1orf111)

Product Code	CSB-MP716346HU
Abbreviation	C1orf111
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T0L3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MENFSLLSIS GPPISSALS AFPDIMFSRA TSLPDIAKTA VPTEASSPAQ ALPPQYQSII VRQGIQNTAL SPDCSLGDTQ HGEKLRNCT IYRPWFSPYS YFVCADKESQ LEAYDFPEVQ QDEGKWDNCL SEDMAENICS SSSSPENTCP REATKKS RHG LDSITSQDIL MASRWHPAQQ NGYKCVACCR MYPTLDFLKS HIKRGFREGF SCKVYYRKLK ALWSKEQKAR LGDRLSSGSC QAFNSPAEHL RQIGGEAYLC L
Source	Mammalian cell
Target Names	SPATA46
Protein Names	Recommended name: Uncharacterized protein C1orf111
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.