



Recombinant Human ATPase family AAA domain-containing protein 3C (ATAD3C)

Product Code	CSB-EP716349HU-B
Abbreviation	ATAD3C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T2N8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSKDALNLAQ MQEQLQLLEQ QSKLKQLVNE DLRKQEEVSQ KHHQTFLESI RAAGTLFGEG FRAFVTDRDK VTATVAGLTL LAVGVYSASN ATAVTGRIYIE ARLGKPSLVR ETSRITVLEA LRHPIQQVSR RLLSRPQDVL EGVVLSPSLE ARVRDIAIMT RNIKKNRGLY RHILLYGPPG TGKTLFAKKL ALHSGMDYAI MTGGDVAPMG REGVTAMHKL FDWANTSRRG LLLFVDEADA FLRKRATEKI SEDLRATLNA FLYRTGQHSN KFMLILASCH PEQFDWAINA CIDVMVHFDL PGQEERARLV RMYLNEYVLK PATEGKRRLK LAQFDYGRKC LEIARLTEGM SCRKIAQLAV SWQATAYASK DGVLTEAMMD ACVQDFVQQH QQMMRWLKG E RPGPEDEQPS S
Source	E.coli
Target Names	ATAD3C
Protein Names	Recommended name: ATPase family AAA domain-containing protein 3C
Expression Region	1-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.