



Recombinant Human Protein lin-9 homolog (LIN9)

Product Code	CSB-MP716370HU
Abbreviation	LIN9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5TKA1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AELDQLPDE SSSAKALVSL KEGSLNNTWN EKYSLSLQKTP VWKGRNTSSA VEMPFRNSKR SRLFSEDDDR QINTRSPKRN QRVAMVPQKF TATMSTPDKK ASQKIGFRLR NLLKLPKAHK WCIYEFYNSN IDKPLFEGDN DFCVCLKESF PNLKTRKLTR VEWGKIRRLM GKPRRCSSAF FEEERSALKQ KRQKIRLLQQ RKVADVSQFK DLPDEIPLPL VIGTKVTARL RGVHDGLFTG QIDAVDTLNA TYRVTFDRTG LGHTHTIPDYE VLSNEPHETM PIAAFGQKQR PSRFFMTPPR LHYTPPLQSP IIDNDPLLGG SPWRSKISGS DTETLGGFPV EFLIQVTRLS KILMIKKEHI KKLREMNTA EKLKSYSMPI SIEFQRRYAT IVLELEQLNK DLNKVLHKVQ QYCYELAPDQ GLQPADQPTD MRRRCESSEEAQ EIVRHANSST GQPCVENENL TDLISRLTAI LLQIKCLAEG GDLNSFEFKS LTDSLNDIKS TIDASNISCF QNNVEIHVAH IQSGLSQMGN LHAFANNTN RD
Source	Mammalian cell
Target Names	LIN9
Protein Names	Recommended name: Protein lin-9 homolog Short name= HuLin-9 Short name= hLin-9 Alternative name(s): Beta subunit-associated regulator of apoptosis TUDOR gene similar protein Type I interferon receptor beta chain-associated pr
Expression Region	2-542
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.