



Recombinant Rat Fin bud initiation factor homolog (Fibin)

Product Code	CSB-EP716390RA-B
Abbreviation	Fibin
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U2T4
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	YF DGPLYPEMSN GTLHHYFVPD GDYEENDDPE KCQLLFRVSD RRRCSQGEGG QASSLLSLTL REEFTVLGRQ VEDAGRVLEG ISKSISYDLD GEESYGKYLR RESHQIGDAY SNSDKSLTEL ESKFKQGQEQ DSRQESRLNE DFLGMLVHTR SLLKETLGIS VGLRDKYELL AHTIRSHGTR LGRLKSDYLE GGGQKMG
Source	E.coli
Target Names	Fibin
Protein Names	Recommended name: Fin bud initiation factor homolog
Expression Region	19-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.