



Recombinant Mouse THO complex subunit 6 homolog (Thoc6)

Product Code	CSB-BP716400MO
Abbreviation	Thoc6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U4D9
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MEHAAPLAVP LGQAEVFAQAL QRLHMTIFSQ SVSPCGKFLA AGNNYGQIAI FSLSAALSSE AKEESKKPVV VFHAHDGPVY SMVSTDRHLL SAGDGEVKGW LWAEILKKGK KELWRRQPPY RTSLEVPEIN ALLLVPKENS LILAGGDCQL HSMDETGFAL TRALRGHTDY IHCLALRERS PEVLSSGGEDG AVRLWDLRIA KEVQTIEVYK HEECSRPHNG RWIGCLATDS DWMVCGGGPA LTLWHLRSST PTTVFPIRAP QKHVTFYQDL ILSAGQGCCV NHWQLSGELK AQVPGSSPGL LSLSLNQQPA APECKVLTAS GNSCRVDVFT NLGYRAFSLF
Source	Baculovirus
Target Names	Thoc6
Protein Names	Recommended name: THO complex subunit 6 homolog Alternative name(s): WD repeat-containing protein 58
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.