



Recombinant Acanthamoeba polyphaga mimivirus Uncharacterized protein R699 (MIMI_R699)

Product Code	CSB-EP716418ADAZ-B
Abbreviation	MIMI_R699
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UNV6
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MEQSNDDNL LVLGIGISVH KTDGVLRF EK YCQAHNLQYM IVGEGKKWNG GNLESEAGGG QKINELLIAL ESIKDNKLIV VCDTYDLIPL SGPEEILRKY RFLTPDNKVV FSSELYCWPD ASLVERYPKV DTKYKYLNSG AFMGYRDDIY EMIKNGVKDR DDDQLFFSIK FIETDKIVLD YKCELFQAMY RCNSDLVVHK NRIFNGYTNS YPVFAHGNP AKKLLNHMEG YFMTEPIDGS SNTINTFKLD NEPKVFFALY VDSNDLSALK QFLGKVASIQ YGNKVIYLYD RSDNEQNRKL IQISYPNYHT GVTKYVFDDF KKSDAQFYFL LEQNCIITKK DILHELMQV KDNHRVISPM IGYEQNSTRT NFWGDIEDGY YKRSENYLDL AKHKVRGLWN VPYVYGVILM HESVVRNWDL SMVKYNDKDM DLCFSLRKHT IFMYMINNNN YGYMV
Source	E.coli
Target Names	MIMI_R699
Protein Names	Recommended name: Uncharacterized protein R699
Expression Region	1-455
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.