



Recombinant Acanthamoeba polyphaga mimivirus Uncharacterized Bro-N domain-containing protein L2 (MIMI_L2)

Product Code	CSB-EP716436ADAZ
Abbreviation	MIMI_L2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UP77
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MRKNLLRILK EKDMSMTKKG YISLYDFVEK IIGSKNPVAY ASKIKDHNPV IINEKEYISP ESCLNILKNA KFAKCRELYN NLQVKDGDKT SIIDVNNNIF QFEGKRFTSF FVDKRDGKWD VWIYGAEVAR FLGYNDCKKA ISIHVESCNR LIFEEIRNNF PIESNSIPKT LDKKTKFINL SGFCNLIHHS KKPFAKIKK WLDDEVIPAL IMDGVYSMQP KELKVKFFYE DNYISDFLLR VIKNRE
Source	E.coli
Target Names	MIMI_L2
Protein Names	Recommended name: Uncharacterized Bro-N domain-containing protein L2
Expression Region	1-246
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.