



Recombinant Acanthamoeba polyphaga mimivirus

Uncharacterized protein L274 (MIMI_L274)

Product Code	CSB-EP716462ADAZ-B
Abbreviation	MIMI_L274
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UPU8
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	SAPCPKKLAD TTCPREPYSF EGMKRQEPIK GRNVGWFQCQ ATMTVHDPVT TAVLVTQPTF TFNLSTLYHN SCGCWKT VTT IPGDAVKMEY GNGLPGQAMY GSCFYATWAD LAGNPVETVG FARNVVG VHR DLIYVSRATN ELVVHQSIYP LGRDQVGDEY VLVRFHDPVT KGIDWVQAVY CTKIQDQLEP LP
Source	E.coli
Target Names	MIMI_L274
Protein Names	Recommended name: Uncharacterized protein L274
Expression Region	21-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.