



Recombinant Acanthamoeba polyphaga mimivirus

Putative ankyrin repeat protein R190 (MIMI_R190)

Product Code	CSB-MP716472ADAZ
Abbreviation	MIMI_R190
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5UQ12
Product Type	Recombinant Protein
Immunogen Species	Acanthamoeba polyphaga mimivirus (APMV)
Purity	≥85% (SDS-PAGE)
Sequence	MEIYADQKYT KLVSHLYYQN SHEIVHSSKI YQLALTNPDI LSQYINYPDE SNESHESHDC HESQNSTHYK FRLMEIYHKM IIDNNMDLL DLLFEIMSDD FRCDDCVILD LAIRHQRTDV FNKYIILGFD LNRLINSSTI LDEIILWYHQ IKKYDTGIEL CDYLIDNGAS ISIHNYCTII NAFDTLNDKY VSFFLNKISL NDLGNLLFWY LRNYYNVNID IVETILSNGI DINNFHEYSY MLIKGFNVPI MNLFIYGLI IHDDVIDDAC RYGNHLLVDY LMEIGHKPSK QIITNVIEH NVNIIKLFVK YNIDLSDIKP PTSPEIIQLV KSLESNGLSV DYVYGYILDK FNLSGKLCI
Source	Mammalian cell
Target Names	MIMI_R190
Protein Names	Recommended name: Putative ankyrin repeat protein R190
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.